

NEW!! **MODEL HS-03 AA**

HEMOSEALER model HS-03 AA is the high frequency tube sealer designed to seal blood bag tubing easily quickly & safely.



SAFE AND SECURE SEALING

The microprocessor in the Power Unit checks the tubing surface in order to avoid sparks leading to blood leakage. The safe sealing procedure makes 3mm wide seal with perforation, which makes it easier to pull the tubes apart.

HEAVY DUTY MODEL

With its powerful high frequency generator, HS-03 AA is suitable not only for routine procedures at donation rooms, but also for the repeated operations at processing rooms without overheating.

BENCH UNIT AND HANDLE UNIT

The Power Unit of HS-03 AA can accept either the Bench Unit, which enables touch-actuated automatic sealing, or the Handle Unit, which has the excellent mobility to reach the sealing part. Both units are with 1.9-meter cable. The sealing electrodes are safely housed and easily accessible for routine cleanings.

ERGONOMIC SEALING HANDLE

Ergonomic Handle with Ceramic Sealing Electrode which minimized the damage of blood cell around sealing portion. Design for easy carry hanging on your shoulder or in your pocket.



Hemosealer MODEL HS-03 AA

Specification

Power Required	110/115/230/240 VAC. 50/60Hz
Battery	Not Available
Power Consumption	Max. 220VA
Continuous operating time on batteries	NIL
Working Frequency	40.68MHz
Sealing Time	Adjustable 0.5 to 5 seconds
Width of Seal	3mm
Dimensions	Power Unit: 290 (L) x 220 (W) x 88 (H) mm Bench Unit: 150 (L) x 60 (W) x 73 (H) mm Handle Unit: 206 (L) x 27 (W) x 35 (H) mm Ergonomic Sealing Handle: 180 (L) x 35 (W) x 145 (H) mm
Weight	Power Unit: 4.1 kgs Bench Unit: 0.9 kgs Handle Unit: 0.3 kgs Ergonomic Sealing Handle Unit: 0.22kgs
Cable Length Ergonomic Sealing Handle Cable Length	1.9m Ca 0.8m spiral (2.3m drawn out)
Electric Shock protective type	Class II type B
Compliance standards	EN60601-1, EN60601-1-2
Code No	Power Unit: EQBHS0101 Bench Unit: EQBHS0102 Handle Unit: EQBHS0103 Ergonomic Sealing Handle Unit: EQBHS0106



The Bridge to Good Health Care

JMS Singapore Pte Ltd

440 Ang Mo Kio Industrial Park 1, Singapore 569620
Tel : +65-6457-1144 | Fax: +65-6459-9564

All copyright reserved