

CLOSED SYSTEM TRANSFER DEVICE

NEO SHIELD

BASIC SYSTEM FOR DRUG COMPOUNDING



**Safe, Reliable
Performance and
Easy to Use**



**Outstanding
Septum Closure**

Adjustable Spike for Various Vials
(13-20mm & 32mm)


Simple Design Supports
Efficient Workflow

SAFE, EASY and EFFICIENT Way of Preparing Antineoplastic Drugs

Vial Spike Designed for Multi-size Drug Vials



■ Procedure (Compounding of Drug with Diluent)



■ Product Line-up



Vial Spike
Code: SV-NSVS
Packing: 200 ea per carton

Vial Spike (32mm)
Code:SV-NSVS32
Packing: 250 ea per carton


Lever Lock
Code: SV-NSLL
Packing: 200 ea per carton

Bag Adaptor
Code: SV-NSBA
Packing: 250 ea per carton

Bag Adaptor - ISO
Code: SV-NSBAI
Packing: 250 ea per carton

Key Point

**Closed
Mechanism of
Lever Lock**



The tip of connection is protected from exposure to environment by the sealing valve. The closed system is maintained during connection and after disengagement.

NEOSHIELD is a registered trademark of JMS Co., LTD.



JMS Singapore Pte Ltd
440, Ang Mo Kio Industrial Park 1
Singapore, 569620
<http://www.jmss.com.sg>

